

LULA Closed Pod System

This system includes a rechargeable device and pre-filled pods, designed to deliver consistent nicotine doses. With child safety features, visible liquid levels, and easy USB-C charging, it helps support your quit journey with simplicity and control.

LULA Therapeutic Vaporiser Device

A reusable, TGA-regulated device with safety features and easy pod swapping.

The device is draw-activated and auto shuts off after 10 minutes. You can find full instructions in our user manual.

Child Lock: Requires 5 quick button presses to activate



LULA Closed Pods

Pre-filled pods that click into your LULA device. Available in 0, 10, 20, 30, 50mg/mL and four flavours with a clear design so you can track usage at a glance.



What We Know About Therapeutic Vaping

Therapeutic vaping is backed by strong evidence and strict regulation. A major review of 90 studies found it can be **twice as effective** as other quit smoking methods, with no serious side effects reported. Therapeutic Vaping Products are regulated by the TGA, so you know exactly what's in them and how much nicotine you're getting. This means better safety, more control and greater confidence on your quit journey. A switch to therapeutic vaping can also save you money, see here:

NVP Cost to Patient*

1 pod per day	\$11
30 pods per month (15 packs)	\$330
24 Week Cessation Program	\$1,980

Cigarette Cost to Patient**

1 pack conventional cigarettes	\$40
One pack a day per month	\$1,200
24 weeks of cigarettes	\$7,200

Savings per month: \$870 | Savings per 24 weeks: \$5,220

*Approximate cost to patient is based on RRP of Closed Pod 20mg/mL. **Calculations are based on avg price of 20 pack cigarettes

How It Works

Insert the Pod

Push pod firmly into the top of the device



Turn On

Press button 5 times quickly



Use the Device

Inhale through the mouthpiece to activate vapour



Turn Off

Press button 5 times or wait 10 minutes for auto shut-off



Using Your LULA Device

Once your device is turned on, it activates automatically when you take a puff through the mouthpiece. You'll feel a gentle vapour as you inhale.

1 Take a steady puff for 2 to 3 seconds, drawing the vapour into your mouth.

2 Remove the device and inhale gently and deeply, then exhale slowly.

Wait at least 20 seconds before taking another puff. This helps you gauge whether your cravings have been satisfied.

Use your LULA Closed Pod System when you feel the urge to smoke. To make each pod last, try to limit your use throughout the day. Each pod delivers around 500 puffs when used in two-second draws. Gradually wean until nicotine free.

Important Information

- Do not leave the device charging unattended
- Keep out of reach of children
- Store away from liquids and direct sunlight
- Dispose of batteries and pods responsibly

Need support to quit?

Visit quit.org.au or download the **My QuitBuddy** app for tips and encouragement.

Questions?

Contact Bay Pharma
1300 938 651 or info@baypharma.com.au - we're here to help.